# ***Video Graphic Printer***

---

## **Instructions for Use** Page 2

Before operating this unit, please read this manual thoroughly and retain it for future reference.

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***UP-880/890CE/890MD***

# Owner's Record

The model and serial numbers are located at the rear. Record these number in the space provided below. Refer to these numbers whenever you call upon your Sony dealer regarding this product.

Model No. \_\_\_\_\_  
Serial No. \_\_\_\_\_

## WARNING

To prevent fire or shock hazard, do not expose the unit to rain or moisture.

To avoid electrical shock, do not open the cabinet. Refer servicing to qualified personnel only.

### For UP-890CE/890MD

#### Symbols on the products

This symbol indicates type B equipment classified in accordance with IEC Publication 601-1 Safety of medical electrical equipment.

This symbol indicates the equipotential terminal which brings the various parts of a system to the same potential.

### For UP-890MD

This symbol is intended to alert the user to the presence of uninsulated "dangerous voltage" within the product's enclosure that may be of sufficient magnitude to constitute a risk of electric shock to persons.

This symbol is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the appliance.

### For the customers in the U.S.A.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

You are cautioned that any changes or modifications not expressly approved in this manual could void your authority to operate this equipment.


The shielded interface cable recommended in this manual must be used with this equipment in order to comply with the limits for a digital device pursuant to Subpart B of Part 15 of FCC Rules.

### For the customers in Europe

#### Important safeguards notices for use in the medical environments

1. All the equipments connected to this unit shall be certified according to Standard IEC601-1, IEC950, IEC65 or other IEC/ISO Standards applicable to the equipments.
2. When this unit is used together with other equipment in the patient area\*, the equipment shall be either powered by an isolation transformer or connected via an additional protective earth terminal to system ground unless it is certified according to Standard IEC601-1.

\* Patient Area

3. The leakage current could increase when connected to other equipment.

### For the customers in the United Kingdom

#### WARNING

#### THIS APPARATUS MUST BE EARTHED

#### IMPORTANT

The wires in this mains lead are coloured in accordance with the following code:

Green-and-yellow: Earth  
Blue: Neutral  
Brown: Live

As the colours of the wires in the mains lead of this apparatus may not correspond with the coloured markings identifying the terminals in your plug proceed as follows:


The wire which is coloured green-and-yellow must be connected to the terminal in the plug which is marked by the letter E or by the safety earth symbol  $\perp$  or coloured green or green-and-yellow.

The wire which is coloured blue must be connected to the terminal which is marked with the letter N or coloured black.

The wire which is coloured brown must be connected to the terminal which is marked with the letter L or coloured red.

# Table of Contents

<b>Introduction</b>	<b>Overview</b>	<b>4</b>
<b>Preparation</b>	<b>Connection</b>	<b>5</b>
	<b>Before Printing</b>	<b>6</b>
	Setting the Type of Paper	6
	Setting the DIP Switches	7
	<b>Loading Paper</b>	<b>10</b>
<b>Operation</b>	<b>Printing</b>	<b>12</b>
	Making Print-outs	12
	Selecting the Printing Direction	13
	Selecting the Printing Size	13
	Adjusting the Contrast and Brightness	15
<b>Others</b>	<b>Precautions</b>	<b>16</b>
	Maintenance	17
	<b>On the Type of Paper</b>	<b>19</b>
	<b>Specifications</b>	<b>20</b>
	<b>Troubleshooting</b>	<b>21</b>
	<b>Location and Function of Parts</b>	<b>22</b>
	Front	22
	Back	23

The UP-880/890CE/890MD is a black and white video graphic printer that can be used to print images displayed on video monitor.

### **Clear, consistent print quality**

- High definition, 10.2 dots/mm printing using a thermal head with high-speed drive IC.
- 256 gradations of black and white.

### **Fast printing**

- You can make a single print-out in about 3.9 seconds in STD and NORM mode.
- You can make a maximum of 11 copies of the same image continuously.

### **Two way printing direction and five printing sizes selectable**

- The printing direction selector on the front panel enables you to print in vertical and horizontal directions.
- The printing size control on the front panel enables you to select five kinds of printing size.

### **DIP switches to optimize the printer**

- You can make print-outs starting either from the bottom or top of the image by setting the DIRECTION DIP switch.
- You can set the print-out aspect ratio to 4:3 or 1:1 by setting the ASPECT DIP switch.
- You can set the range to be printed by setting the SCAN DIP switch.
- You can save your paper by setting the POSTFEED DIP switch (paper saving function).

### **Automatic video signal discrimination**

The type of input signal, black and white (EIA or CCIR) or input color (NTSC or PAL), is automatically discriminated and printed in same duration and size.

### **Alarm buzzer**

The alarm buzzer prevents you from making any misoperation.


### **Easy and quick paper loading**

You can load paper just by opening the paper lid with the OPEN/CLOSE button and placing the paper roll.

# Connection

## Notes

- Turn off the power to each device before making connection.
- Connect the AC power cord last.

\* The RM-91 remote commander is supplied with UP-890MD only.


# Before Printing

## Setting the Type of Paper

Set the PAPER TYPE selector to the type of paper to be used.

The use of paper other than Sony may result in reduced printer performance and poor print quality.

Type of paper	PAPER TYPE switch position
UPP-110S	I (Normal)
UPP-110HD	II (High density)
UPP-110HA	IV (Enhanced)

### When you use the UPP-110HA or UPP-110HD

When you set the PAPER TYPE selector to II or IV, set the density gradation with the GAMMA selector.


I: Soft gradation

II: Standard

IV: Hard gradation

## Setting the DIP Switches

Set the DIP switches according to the required print mode. Before setting the DIP switches, turn the power off. Change the settings using a small pointed tool such as a small screwdriver. The factory settings are as follows.

### 1 INTERRUPT ON/OFF switch

To interrupt the printing under way and print a new picture when you press the PRINT button during printing, set this switch to ON.

To disregard that the PRINT button is pressed during printing and continue the printing under way, set to OFF.

If you press the PRINT button during printing in OFF mode, the alarm buzzer will sound.

### 2 POSTFEED ON/OFF switch


To feed out extra blank paper once a picture has been printed, set this switch to ON.

To save paper by feeding only a short length of paper after printing a picture, set to OFF. You can make more print-outs per roll of printing paper, but you have to take out and cut the paper yourself.

### 3 ASPECT 4:3/1:1 switch

Normally keep this switch set to 4:3. When the aspect ratio of the video signal is 1:1, set to 1:1.

The print-out will be longer than a print-out made at 4:3.

### 4 MEMORY FRAME/FIELD switch

Normally keep this switch set to FRAME (ON). When printing fast-moving pictures (such as a ball being thrown), the print-out may blur. If this happens, set to FIELD. The print-out definition will be poorer but less blurred.

### 5 IMAGE POSI/NEGA switch

Normally keep this switch set to POSI (ON). To make negative print-outs, set to NEGA (OFF).


### 6 MIRROR NORM/REV switch

Normally keep this switch set to NORM (ON). To print the right and left sides reversed, set to REV (OFF).

### 7 DIRECTION NORM/REV switch

Selects whether the top or bottom of the screen is to be printed first.


Normally keep this switch set to NORM (ON). Printing is done from the bottom of the screen. To start printing from the top of the screen, set to REV (OFF).

### **8 9 SCAN NORM/WIDE 1/WIDE 2 switch**

Sets the print-out range. The print-out range is widened in the NORM, WIDE 1, and WIDE 2 order.

To print only the image displayed on the standard screen size of the video monitor, set the SCAN switch **9** to NORM (OFF). To print when the signal scans beyond the edge of the standard monitor screen, set the SCAN switch **9** to WIDE 1 (ON) or the SCAN switch **8** to WIDE 2 (OFF). When you set the SCAN switch **8** to the WIDE 2 position, WIDE 2 is selected regardless of the setting of the SCAN switch **9** position.

### **10 RESERVED switch**

Keep this switch set to ON.

### **11 INPUT B&W/COLOR switch**

Set this switch to B & W (ON) when the signal to be printed is black and white. Set to COLOR (OFF) when the signal is color.

### **12 75Ω ON/OFF switch**

Set this switch to OFF when a video monitor or other video equipment is connected to the VIDEO OUT connector.

Set to ON when nothing is connected to the VIDEO OUT connector.

When you connect two printers to one video equipment, set the 75Ω switch of one of the printer to ON, and the other to OFF.


# Loading Paper

## Notes


- Do not fold the paper or touch the printing surface. Dust on the printing surface will result in poor print quality.
- After loading the paper roll, pull out and cut off the first 15 to 20 cm (6 to 7 7/8 inches) to remove any slack.
- Use only UPP-110 series paper (p. 19).
- Set the PAPER TYPE selector according to the paper type (p. 6).

## Loading


- 1 Press the power ON/OFF switch to turn on the printer.

- 2 Press the OPEN/CLOSE button to open the paper lid.



- 3 Place the paper roll in the printer.


Place the paper with the thermo-sensitive side (printing side) up.

- 4** Pull out the first 15 to 20 cm (6 to 7 7/8 inches) of the paper to remove any slack in the roll.

- 5** Press the OPEN/CLOSE button to close the paper lid. You can also close the paper lid simply by pushing it.

**CAUTION**

Keep fingers clear of paper lid assembly and paper cutting blade when paper lid is closing.


# Printing

## Before making print-outs


- Are the connections correct? (p. 5)
- Is the paper roll loaded properly? (p. 10).
- Is the paper type set correctly? (p. 6).
- Are the DIP switches set correctly? (p. 7 – 9)
- Is the print source being input?

## Making Print-outs

- 1 Press the power ON/OFF switch to turn on the printer.  
The power indicator lights.


- 2 Make sure that the PAPER EMPTY indicator is not lit.  
If lit, load paper.
- 3 Select the printing direction and size.  
See "Selecting the Printing Direction" and "Selecting the Printing Size" on the next page.
- 4 When the picture you want to print is on the video monitor, press the PRINT button.  
The printer makes a print-out of the picture displayed at the instant you press the PRINT button.


## Selecting the printing direction

You can select the vertical or horizontal direction using the STD/SIDE selector.  
 To print in the vertical direction, set to STD.  
 To print in the horizontal direction, set to SIDE.

## Selecting the printing size

You can print in small or large size.  
 Also you can enlarge the center of the picture by 1.5 or 2 times.

Control position (mode)	Printing size	
	Vertical (STD)	Horizontal (SIDE)
NORMAL	ABC	ABC
SMALL	ABC ABC	ABC ABC
LARGE	ABC	ABC
ZOOM 1.5X	ABC	ABC
ZOOM 2X	ABC	ABC

### To print in SMALL mode

Press the PRINT button twice. When you press the PRINT button once, the buzzer sounds. The printer starts printing after the PRINT button is pressed twice.

### Stopping printing midway

Press the OPEN/CLOSE button while printing or while copying. The printer stops printing.

#### To stop printing and print another picture displayed on the video monitor

To do this, the DIP switch **1** (INTERRUPT) must be set to ON (p. 7). Press the PRINT button while printing or copying. The printer stops printing and starts printing the picture displayed at the instant you press the PRINT button.

### Making copies of the last print-out

Press the COPY button. The printer makes a copy of the last print-out. The last print-out is retained in the printer's memory until you press the PRINT button again or turn the power off.

#### To copy in different sizes

You can copy the last print-out in different sizes. Before pressing the COPY button, select the printing size as described in "Selecting the Printing Size".

#### Notes

- If you press the COPY button immediately after turning the power on, the alarm buzzer will sound as nothing is stored in memory.
- In SMALL mode (p. 13), if you press the COPY button after you have pressed the PRINT button only once, the alarm buzzer will sound and the printer will not copy.

#### To make multiple copies of the same print-out

Press the COPY button as many times as necessary (maximum 11 copies including the first print-out) while printing or copying the first print-out. Each time you press the COPY button, the short buzzer sounds.


#### To stop copying midway

Press the OPEN/CLOSE button.

# Adjusting the Contrast and Brightness


You can adjust the contrast and brightness of the print-out.

- 1 Set the THRU/EE selector to EE.  
You can check the adjusted picture on the video monitor.


- 2 Adjust the brightness with the BRIGHT control and the contrast with the CONTR control while watching the picture on the video monitor.

## Adjusting the contrast

## Adjusting the brightness

## To directly input the video signal from the video equipment, which is connected to the printer, to the video monitor

Set the THRU/EE selector to THRU. The video signal is directly input to the video monitor without being processed by the printer's circuitry.

# Precautions

---

## On the safety

- Check the operating voltage before operation.  
Operate the unit only with a power source specified in “Specifications”.
- Stop operation immediately if any liquid or solid object falls into the cabinet.  
Unplug the unit and have it checked by qualified personnel.
- Unplug the unit from a wall outlet if you will not be using it for a long time.  
Disconnect the power cord by grasping the plug. Never pull the cord itself.
- Do not disassemble the cabinet. Refer servicing to qualified personnel only.
- Do not touch the cutting blade of the printer.
- Keep fingers clear of paper lid assembly and paper cutting blade when paper lid is closing.
- Connect the power plug of the printer to a wall outlet having protective earth terminal. The safety earth should be properly established.

---

## On operation

Do not turn the power off while the printer is printing. The thermal head may be damaged.

---

## On printer carriage

Do not carry and move the printer when the paper roll is placed in the printer. Doing so may cause malfunction.

---

## On installation

- Place the printer on a level and stable surface during operation.
- Do not install the printer near heat sources. Avoid locations near radiators or air ducts, or place subject to direct sunlight or excessive dust, humidity, mechanical shock or vibration.
- Provide adequate air circulation to prevent heat build-up. Do not place the printer on surfaces such as rugs, blankets, etc., or near materials such as curtains and draperies.

# Maintenance


## Cleaning the cabinet

Do not use strong solvents to clean the printer. Thinner or abrasive cleansers will damage the cabinet.


## Cleaning the thermal head

If the print-out is dirty or white stripes appear on the print-outs, clean the thermal head using the supplied cleaning sheet.


- 1 Press the power ON/OFF switch to turn on the printer.


- 2 Press the OPEN/CLOSE button to open the paper lid.

- 3** Insert the cleaning sheet, with the black surface facing down, into the groove in the paper lid.

- 4** Press the OPEN/CLOSE button and keep it pressed.  
The paper lid closes and the printer starts cleaning the head.  
When the buzzer sounds and the printer starts ejecting the cleaning sheet, release the OPEN/CLOSE button.

- 5** Remove the cleaning sheet.

### Notes

- Do not press the PRINT or COPY button while the cleaning sheet is in the printer.
- Clean the head only when necessary. If you clean the head too often, it may cause malfunction.

# On the Type of Paper

---

## Type of paper

- Use only the Sony UPP-110 series paper. The use of other paper may result in reduced printer performance and poor print quality.
- The following types of paper are available.

Printing density	Type of paper
TYPE I (Normal)	UPP-110S
TYPE II (High density)	UPP-110HD
TYPE IV (Enhanced)	UPP-110HA

---

## Storing paper

- Store unused or printed paper in a cool, dark place (below 30°C or 86°F). We recommend that you store printed paper in a polypropylene pouch.
- Do not store unused or printed paper in hot or humid place.
- Do not leave unused or printed paper in direct sunlight or other bright place for extended periods.
- Do not allow any volatile organic solvent or vinyl chloride to touch the printed paper. Alcohol, plastic tape or film will fade the print-out.
- To attach printed paper to another piece of paper, use double-sided adhesive tape, or water-based or solid glue.
- Do not stack printed paper on or under a diazo copy sheet. The print-out may become discolor in black.

# Specifications

## Thermal head \*

Thin-film thermal head (with built-in drive IC) 1024-dot drive

## Gradation

256

## Resolution (in WIDE 1 mode)

EIA: 970 x 490 dots

CCIR: 970 x 582 dots

## Print size (in NORM and WIDE 1 mode)

### STD mode

EIA: 95 x 72 mm

CCIR: 95 x 71 mm

### SIDE mode

EIA: 127 x 96 mm

CCIR: 127 x 95 mm

## Printing speed (in STD and NORM mode)

About 3.9 seconds/screen (aspect ratio 4:3)

## Picture memory

768 K x 6 bits

## Input/output connectors

### VIDEO IN (BNC)

EIA or CCIR

Composite video signals

1.0 V<sub>p-p</sub>, 75 ohms/high-impedance (EIA/

CCIR automatically discriminated)

### VIDEO OUT (BNC)


EIA or CCIR

Composite video signals

1.0 V<sub>p-p</sub>, 75 ohms, loop-through/EE

switchable

### REMOTE (stereo minijack)

1 GND

2 PRINT SIGNAL (TTL)

Input of LOW pulse over 100 msec.  
initiates print.

3 PRINT BUSY (TTL)

Goes HIGH during printing.

## Power requirements and consumption

UP-880: 120 V AC, 50/60 Hz, 1.8 A

UP-890CE/890MD:

120 V AC, 50/60 Hz, 1.5 A

220 to 240 V AC, 50/60 Hz, 0.8 A

## Dimensions

Approx. 154 x 106 x 303 mm (w/h/d) (6 1/8 x 4 1/4 x 12 inches)

## Mass

Approx. 3.5 kg (7 lb 11 oz), Main unit only

## Supplied accessories

Paper roll (UPP-110HA) (1)

BNC – BNC connecting cable (1)

AC power cord (1)


Head cleaning sheet (1)

Remote commander RM-91 (1) supplied  
with UP-890MD only

Design and specifications are subject to change  
without notice.

# Troubleshooting


Symptom	Cause/remedy
White specks on first few print-outs.	When printing with a newly inserted roll of paper, dust on the surface of the paper may cause white specks on the print-outs. → Feed the paper by pressing the OPEN/CLOSE button until clean paper appears.
Printing does not start when you press the PRINT button.	<ul style="list-style-type: none"> <li>Paper is not fed. → Is the paper slack? → Is the power turned on? → Are all connections correct? (p. 5) → Did you press the PRINT button twice in SMALL mode?</li> <li>When the alarm buzzer sounds: → Has the thermal head overheated? → Is the video signal of the picture input? → Is the paper loaded correctly?</li> <li>Paper is fed, but printing does not start. → Is the paper loaded with the thermo-sensitive side up?</li> </ul>
Black borders or missing portions around the print-out.	This may result according to the video signal input to the printer. → Change the setting of the SCAN switches (DIP switches [8], [9]). (p. 9)
Paper jam	<ul style="list-style-type: none"> <li>Open the paper lid by pressing the OPEN/CLOSE button, then pull the jammed paper slowly and remove it.</li> <li>There is condensation within the unit. → Moving the unit suddenly from a cold place to a warm place often results in condensation forming. In the event of condensation forming, remove the paper, turn off the power and leave the unit for about one to two hours.</li> </ul>
Print-out is dirty.	The thermal head is dirty. → Clean the thermal head with the supplied head cleaning sheet. (p. 17)
The printer stops printing when it prints continuously black pictures.	This is likely to occur when the printer prints continuously 15 or more black pictures. In such a case, the buzzer sounds. This is because that the protective circuit works against heat build-up of the thermal head. Stop printing for a while.
White lines or small letters on the screen are not printed clearly.	Is the INPUT switch (DIP SWITCH [11]) set to B & W when the input signal is a black and white signal? (p. 9)
Small squares appear over the whole screen.	Is the INPUT switch (DIP switch [11]) set to COLOR when the input signal is a color signal? (p. 9)
The print-out is too dark or too light.	<ul style="list-style-type: none"> <li>Is the 75Ω switch (DIP switch [12]) set correctly? (p. 9)</li> <li>Is the GAMMA selector set correctly? (p. 6)</li> </ul>
The print-out seems stretched.	The ASPECT switch (DIP switch [3]) is set to 1:1. → Set to 4:3. (p. 8)

# Location and Function of Parts

For details, refer to the pages indicated in parentheses.

## Front

**1 Power ON/OFF switch and indicator**  
Turns the power on. The indicator is lit while the power is on.

**2 Printing size selector (13)**  
Selects the printing size.

**3 STD (standard)/SIDE selector (13)**  
Selects the printing direction.

**4 CONTR (contrast) control (15)**  
Adjusts the contrast of the print-outs.

**5 THRU/EE selector (15)**  
Selects the video signal output from the VIDEO OUT connector.  
**THRU:** Input signals are directly output to the video monitor.  
**EE:** Input signals are output to the video monitor after being processed by the printer's circuitry.

**6 BRIGHT (brightness) control (15)**  
Adjusts the brightness of the print-outs.

**7 COPY button (14)**  
Prints another copy of the previous print-out.


**8 PRINT button (12)**  
Prints the picture currently displayed on the video monitor. The picture displayed when you press the PRINT button is stored in memory.

**9 PAPER EMPTY indicator**  
Lights when the printer is out of paper.

**10 OPEN/CLOSE button (10, 14)**  
Opens or closes the door. Also, stops printing midway.

**11 Paper feeder and cutter**  
Cuts the printing paper.

## Back

### **1 REMOTE connector (5)**

Connect the RM-91 remote commander for controlling print operation from a distance. The RM-91 remote commander is supplied with UP-890MD only.

### **2 GAMMA selector (6)**

Changes the print mode to that for high-density printing paper. The selector is effective when the PAPER TYPE selector is set to II or IV.

### **3 PAPER TYPE selector (6)**

Sets the type of paper.

### **4 DIP SW (switches) (7 – 9)**

Sets the print modes and functions.

### **5 VIDEO IN (input) connector (BNC type) (5)**

Connect to the video output connector of the video equipment.

### **6 VIDEO OUT (output) connector (BNC type) (5)**

Connect to the video input connector of the video monitor. The output signal type depends on the setting of the THRU/EE selector.

### **7 AC IN (AC power input) connector (5)**

Connect to a wall outlet using the supplied AC power cord.

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